

The background features a large, light gray hexagon. Inside it is a smaller, light green hexagon. In the center of the green hexagon is a circular warning sign. The sign has a green border and a green background with a white exclamation mark. The text "Calculating the limit of detection (LOD) for a specific analyte with a spectrometer" is overlaid on the center of the image.

Calculating the limit of detection (LOD) for a specific analyte with a spectrometer

Created by: Francheska Colon, Naihomy Tirado, Edwin Caballero

Manager: Samuel Hernandez-Rivera

SOP-01	Francheska Colon Naihomy Tirado, Edwin Caballero		University of Puerto Rico at Mayagüez
Effectivity: July/10/2022	Calculating the limit of detection (LOD) for a specific analyte with a spectrometer		Revised by:
Revised:			Approved by:

This SOP uses the following:

- Equipment: Desired spectrometer.
- Instruments: Volumetric flask, pipettes, and sample container.
- Programs: Microsoft Excel.
- Reagents: Analyte and solvent of interest.

I. PLANNING INITIAL CONCENTRATION RANGE SCREENING

1. **Research** the LOD of the analyte of interest for the phenomenon and technique desired.

E.g., LOD for analyte “X” is 0.400 M.

2. **Determine** the upper limit (UL) by multiplying the LOD by 5.

$$E.g., UL = LOD \cdot 5 \rightarrow UL = 0.400 M \cdot 5 \rightarrow UL = 2.00 M$$

3. **Determine** the increments by subtracting the UL from the lower limit (LL) and dividing by 4.

$$Increments = \frac{UL - LL}{4} \rightarrow Increments = \frac{(2.00 - 0.400) M}{4} \rightarrow Increments = 0.400 M$$

2.00 M	1.60 M	1.20 M	0.800 M	0.400 M
-------------------	-------------------	-------------------	--------------------	--------------------


4. **(OPTIONAL)** IF LOD of the analyte of interest is not found, use the range from 50% to 10% at 10% increments.

50%	40%	30%	20%	10%
------------	------------	------------	------------	------------


5. **Determine** the volume of the available flasks and the amount per volume.




6. **Determine** the final volume of your solution (V_2), i.e., the flask that will be used for the solution.



7. **Find** the initial concentration (C_1) of the analyte of interest.



8. **Determine** the desired concentration of your solution (C_2).



9. **Determine** the volume needed for analyte with concentration C_1 to reach the desired concentration (V_1) with a final volume of V_2 .

$$V_1 = \frac{C_2 \cdot V_2}{C_1}$$

10. (OPTIONAL IF AN INTEGER V1 VALUE IS NEEDED) Determine the final concentration C2 for the analyte with an initial concentration C1 after adding a specific initial volume V1, close to the value obtained from step 6, and completing to a final volume of V2.

$$C_2 = \frac{C_1 \cdot V_1}{V_2}$$

11. Repeat steps 3-7 for the other solutions.


IF an integer value of V1 cannot be obtained for a solution with the initial C1 concentration of the analyte of interest, use the final concentration C2 of a solution as the initial concentration C1 to determine the volume.

Solution used to create other solutions might be consumed, create more, or change V2 to a bigger value.

E.g.

Solution	C1 [%v/v]	C2 [%v/v]	V2 [mL]	V1 [mL]
blank	0%	0%	10.00	10.00
1	99.00%	31.680%	25.00	8.00
2	31.68%	25.344%	10.00	8.00
3	31.68%	22.176%	10.00	7.00
4	99.00%	19.800%	50.00	10.00
5	19.80%	15.840%	25.00	20.00
6	15.84%	12.672%	10.00	8.00
7	99.00%	9.900%	10.00	1.00


12. Prepare solutions by adding the initial volume of analyte (V1) of C1 concentration into the volumetric flask with the final volume (V2) to obtain the final concentration C2.




II. SCREENING LOD CONCENTRATION RANGE FOR ANALYTE

1. **Acquire** spectra for each concentration and blank.

10 measurements per sample. Average the measurements.



2. **Determine** the peak with intensity proportional to the analyte concentration.



3. **Limit** spectral data to the frequencies of the chosen spectral band.



frequency	958	962	966	970	974	978	982	986	990	994
intensity	0.2666	0.331	0.396	0.45	0.469	0.436	0.337	0.291	0.257	0.22

4. **Calculate** the area of the chosen peak with the first frequency (ν_1) and intensity (I_1) as well as the second frequency (ν_2) and intensity (I_2) with the following equation.

$$A = (\nu_2 - \nu_1) \cdot \left(\frac{I_1 + I_2}{2} \right)$$

frequency	958	962	966	970	974	978	982	986	990	994
intensity	0.2666	0.331	0.396	0.45	0.469	0.436	0.337	0.291	0.257	0.22

Diagram illustrating the calculation of area for a peak. Red arrows point from the frequency values ν_1 and ν_2 to the corresponding intensity values I_1 and I_2 in the table. The values 958, 962, 0.2666, and 0.331 are circled in red.

5. **Repeat** step 4 until the area for each frequency is calculated (first frequency does not contain area).

frequency	958	962	966	970	974	978	982	986	990	994
area		1.195	1.455	1.693	1.838	1.81	1.545	1.255	1.095	0.954

Diagram illustrating the calculation of area for a peak. Red arrows point from the frequency values ν_1 and ν_2 to the corresponding area values A_1 and A_2 in the table. The values 958, 962, 1.195, and 1.455 are circled in red.

6. **Add** the area for each frequency to obtain the area of the peak (PA).



$$PA = \sum_{i=1}^n [A_i]$$

$$\rightarrow PA = 0 + 1.195 + 1.455 + 1.693 + 1.838 + 1.810 + 1.545 + 1.255 + 1.095 + 0.954$$

$$\rightarrow PA = 12.841$$


7. **Repeat** steps 4-6 for all samples.

8. Plot peak area vs concentration.



9. Calculate slope and intercept error by using the = LINEST function.


Select the 'y' values, then the 'x' values, write "TRUE" twice, and click enter. A range of values will appear below (2 columns and 5 rows).



Names will not appear.

Slope	0.1068	19.458	Intercept (b)
Standard Error m	0.1184	6.832	Standard Error b
R2	0.092	10.7868	Standard Error Reg
Fisher's F	0.81	8	dof
SS Reg	94.8	930.8	SS Res


10. Add trendline by selecting the data points, right click, and choose “Add Trendline...”.



11. Write linear equation with R^2 and errors.

$$y = \text{slope}(\pm \text{standard error } m) + \text{intercept}(\pm \text{standard error } b)$$


$$R^2$$




12. Calculate standard deviation (s) for each concentration.

$$= STDEV(\text{Range of values})$$




13. Add error bars by selecting the data points, selecting the “Chart Design” tab, “Add Chart Element” button, “Error Bars” list element, and “More Error Bars Options” option.



14. Use standard deviation by going to “Custom” in Error Amount and “Specify Value”.



15. Select standard deviation values in both positive and negative error values.



16. Calculate relative standard deviation (RSD) for the first row by dividing the average peak area with the standard deviation.

$$RSD = \frac{s}{\bar{x}} \times 100\%$$

desired average value

desired standard deviation value

desired RSD value (first row)

Concentration, %v/v	Peak Area 1120 cm ⁻¹	Standard Deviation n	Coefficient of Determination (R ²)	Relative Standard Deviation (RSD)
4.5	16.000	1.5	-	=D5/C5
15.5	18.000	1.6	1.000	8.9%

17. Repeat step 9 for the next peak.

desired average value

desired standard deviation value


desired RSD value (first row)

Concentration, %v/v	Peak Area 1120 cm ⁻¹	Standard Deviation n	Coefficient of Determination (R ²)	Relative Standard Deviation (RSD)
4.5	16.000	1.5	-	9.4%
15.5	18.000	1.6	1.000	=D6/C6


18. Repeat step 10 for row with peak values.

Concentration, %v/v	Peak Area 1120 cm ⁻¹	Standard Deviation n	Coefficient of Determination (R ²)	Relative Standard Deviation (RSD)
4.5	16.000	1.5	-	9.4%
15.5	18.000	1.6	1.000	8.9%
24.5	24.000	1.8	0.890	7.5%
35.5	26.000	2.1	0.935	8.1%
44.5	47.000	7.0	0.785	14.9%
55.5	18.000	2.9	0.197	16.1%
64.5	15.000	4.0	0.018	26.7%
75.5	19.000	4.2	0.000	22.1%
84.5	25.000	5.6	0.003	22.4%
95.5	40.000	12.0	0.092	30.0%



19. Graph relative standard deviation (RSD) vs concentration.



20. Determine the concentration range by analyzing the RSD values.



Values above 10% can be discarded from the desired LOD range. E.g., the concentration range for this data set is 9.9 – 35%v/v.



III. CALCULATING LOD FOR SPECIFIC ANALYTE

1. **Calculate** the increment concentrations for the range determined in the previous section (II).


E.g., calculate the increments for the range 9.9 – 35.5%v/v for 5 samples.

$$\text{Increments} = \frac{UL - LL}{n - 1} \rightarrow \text{Increments} = \frac{(35.5 - 9.9) \%}{5 - 1} \rightarrow \text{Increments} = 6.4 \%$$


9.9 %	16.3 %	22.7 %	29.1 %	35.5 %
----------	-----------	-----------	-----------	-----------

2. **Prepare** solutions as shown in section I.

$$V_1 = \frac{C_2 \cdot V_2}{C_1}$$




3. **Acquire** spectra for each sample.
4. **Repeat** step 3 until 10 measurements are acquired for each concentration and 20 measurements for the blank.





5. **Calculate** the area for the peak of interest for each concentration as in step II.

$$A = (v_2 - v_1) \cdot \left(\frac{I_1 + I_2}{2} \right)$$



6. **Plot** peak area vs concentration.




7. **Register** the slope (m) for the calibration curve.

$$y = 0.0069(\pm 0.0004) + 0.0609(\pm 0.0103) \rightarrow m = 0.0069$$


8. **Calculate** the average for all 20 peak area measurements for the blank sample.

$$\bar{S}_{bl} = \frac{\sum_{i=1}^n [PA_i]}{n}$$



9. **Calculate** the standard deviation for all 20 peak area measurements for the blank sample.

$$s_{bl} = \sqrt{\frac{\sum_{i=1}^n [(PA_i - \overline{PA})^2]}{n - 1}}$$



10. **Calculate** the LOD with the average (\bar{S}_{bl}), standard deviation (s_{bl}), and regression slope (m).

$$LOD = \frac{(\bar{S}_{bl} + 3.3 \cdot s_{bl}) - \bar{S}_{bl}}{m} \rightarrow LOD = \frac{(0.018 + 3.3 \cdot 0.006) - 0.018}{0.0069} \rightarrow LOD = 2.9\%$$

11. **Calculate** the LOQ with the average (\bar{S}_{bl}), standard deviation (s_{bl}), and regression slope (m).

$$LOQ = \frac{(\bar{S}_{bl} + 10 \cdot s_{bl}) - \bar{S}_{bl}}{m} \rightarrow LOQ = \frac{(0.018 + 10 \cdot 0.006) - 0.018}{0.0069} \rightarrow LOQ = 8.7\%$$

Advisor Signature

Co-Advisor Signature
